

Alignment of Prox-seq sequencing data

Hoang Van Phan

Tay Lab, The University of Chicago

Alignment workflow

1. Extract antibody identity: ReadAlignmentDropSeq or ReadAlignmentSmartSeq.
2. Merge cell barcode (if needed): BuildCellBarcodes.
3. Cell barcode correction (if needed): CellBarcodeCorrection.
4. Identify true cell barcodes using the knee plot (if needed): ReadcountHistogram.
5. Merge UMI: UMIMerging.
6. Extract digital count: DigitalCount.

Below is an example usage of the alignment program, and detailed descriptions of each step.

Antibody barcode alignment

For the first step, we have to extract the antibody identities, cell barcodes, and UMIs from the sequencing reads. There are currently two options, depending on whether the Drop-seq or plate-based version of Prox-seq was used.

Drop-seq:

```
java -jar PLA_alignment.jar ReadAlignmentDropSeq \  
R1=read1.fastq.gz \  
R2=read2.fastq.gz \  
ABfile=barcode_cocktail.csv \  
O=ReadAlignment_out.txt.gz \  
SUMMARY=ReadAlignment_summary.txt \  
HEADER=TRUE
```

Arguments:

R1 and R2: paths to the read 1 and read 2 files (fastq.gz format). Read 1 contains the cell barcodes and the UMIs, just like Drop-seq's read 1 files.

ABfile: path to the comma-separated file containing the protein target names (first column) and their barcode sequences (second column).

Example: CD3,AGTCGACA

O: path to store the alignment results (txt.gz format).

SUMMARY (optional): path to store the summary file (txt format). The default is the current working directory.

HEADER (optional): indicates if ABfile contains a header row. The default is FALSE.

Plate-based:

```
java -jar PLA_alignment.jar ReadAlignmentSmartSeq \  
R1List=R1List.csv \  
ABfile=barcode_cocktail.csv \  
O=ReadAlignment_out.txt.gz \  
SUMMARY=ReadAlignment_summary.txt \  

```

HEADER=TRUE

Arguments:

R1List: path to the csv file containing the paths to single cell read 1 files, and their cell id.

The csv file must NOT contain a header row.

Example: /path/to/cell_1_read1.fastq.gz,cell_1

ABfile, O, SUMMARY, and HEADER options: identical to ReadAlignmentDropSeq.

Cell barcode correction

This step is only required for Drop-seq data. The user either supplies a list of “correct” reference cell barcodes (CellBarcodeCorrection), either provided by Drop-seq alignment tools (output of BAMTagHistogram function) or generated using BuildCellBarcodes.

CellBarcodeCorrection

```
java -jar PLA_alignment.jar CellBarcodeCorrection \
I=ReadAlignment_out.txt.gz \
O=cellbarcode_out.txt.gz \
CELL_BC_LIST=referencebarcodes.txt.gz \
READCOUNT_CUTOFF=100 \
SUMMARY=cellbarcode_summary.txt \
HEADER=TRUE
```

Arguments:

I: path to alignment results (output of ReadAlignmentDropSeq).

O: path to store output file with cell barcode corrected (txt.gz format).

CELL_BC_LIST: path to a tab-separated list of cell barcodes (second column) and their read counts (first column). This is the file produced by Drop-seq alignment tools (output of BAMTagHistogram function), or by the function BuildCellBarcodes below.

READCOUNT_CUTOFF (optional): the minimum read count for a reference cell barcode to be used. The default is 100.

SUMMARY (optional): path to store summary file. The default is the current working directory.

HEADER (optional): indicates if CELL_BC_LIST has a header row to be discarded. The default is FALSE.

BuildCellBarcodes

```
java -jar PLA_alignment.jar BuildCellBarcodes \
I=ReadAlignment_out.txt.gz \
O=cellbarcode_out.txt.gz \
SUMMARY=cellbarcode_summary.txt
```

Arguments:

I: path to alignment results (output of ReadAlignmentDropSeq).

O: path to store output file with cell barcode corrected.

SUMMARY (optional): path to store summary file. The default is the current working directory.

This function builds a list of reference cell barcodes by merging all cell barcodes that are 1 Hamming distance apart, and keep the cell barcode with the highest read count. If a cell

barcode matches with more than 2 cell barcodes at 1 Hamming distance, then this barcode is discarded.

NOTE: the output file is tab-separated, and has a header row just like the output of BAMTagHistogram from Drop-seq alignment tools.

Knee plot

This step is only optionally used for Drop-seq data. Here, we export the read count of each single-cell barcode.

```
java -jar PLA_alignment.jar ReadcountHistogram \
I=cellbarcode_out.txt.gz \
O=cellbarcode_readcounts.txt.gz
```

Arguments:

I: path to the cell barcode corrected file.

O: path to export the tab-separated cell barcodes (second column) and their corresponding read counts (first column) (txt.gz format).

Then, we can use the output cellbarcode_readcounts.txt.gz to plot the knee plot to identify the good single-cell barcodes (similar to Drop-seq tools).

UMI merging

Now, we have to merge PLA products coming from the same single cell with similar UMI sequences together. Here, reads with duplicated UMIs, and reads with UMIs that are 1 Hamming distance from another UMI with a higher read count are discarded. This step is performed for both Drop-seq and plate-based versions.

Drop-seq:

```
java -jar PLA_alignment.jar UMIMerging \
I=cellbarcode_out.txt.gz \
O=umi_out.txt.gz \
SUMMARY=umi_summary.txt
```

Plate-based:

```
java -jar PLA_alignment.jar UMIMerging \
I= ReadAlignment_out.txt.gz \
O=umi_out.txt.gz \
SUMMARY=umi_summary.txt
```

Arguments:

I: path to the input file.

O: path to export the UMI merging result (txt.gz format).

SUMMARY (optional): path to store summary file. The default is the current working directory.

Export to digital count matrix

Finally, we export a count matrix containing the UMI counts of all single-cell PLA products. The count matrix format is PLA product by single cells.

Drop-seq:

```
java -jar PLA_alignment.jar DigitalCount \  
I=umi_out.txt.gz \  
O=count_matrix.txt.gz \  
CELL_BC_LIST=knee_barcodes.txt \  
HEADER=FALSE \  
SUMMARY=count_summary.txt \  
DUPLICATE_EXPORT=duplicate_export.txt.gz \  
REMOVE_DUPLICATE=TRUE
```

Plate-based:

```
java -jar PLA_alignment.jar DigitalCount \  
I=umi_out.txt.gz \  
O=count_matrix.txt.gz \  
CELL_BC_LIST=NONE \  
SUMMARY=count_summary.txt \  
DUPLICATE_EXPORT=duplicate_export.txt.gz \  
REMOVE_DUPLICATE=TRUE
```

Arguments:

- I: path to the UMI merging result.
- O: path to store the tab-separated digital count matrix. Rows are PLA products, columns are single-cell barcodes (txt.gz format).
- CELL_BC_LIST (optional): a list of chosen cell barcodes (from knee plot) (txt format). If NONE (default), export all available cell barcodes.
- HEADER (optional): whether the CELL_BC_LIST has a header row. The default is FALSE.
- SUMMARY (optional): path to store the summary file. The default is the current working directory.
- DUPLICATE_EXPORT (optional): path to store the list of duplicated PLA products across single cells. These are the PLA products that have the same UMI in more than one cell. The default is the current working directory.
- REMOVE_DUPLICATE (optional): whether to remove duplicated PLA products across single cells. The default is FALSE.